


Initiated: 2/15/12  
 Revised: 6/1/15  
 Revision: 2

**MILWAUKEE COUNTY EMS  
 PRACTICE GUIDELINE  
 12-LEAD ECG ACQUISITION**

Approved: M. Riccardo Colella, DO, MPH, FACEP  
 Page 1 of 1



- Notes:**
- Leads (electrodes) should not be removed; wires can be removed if necessary but prefer to simply unplug cable from monitor and reconnect in rig so serial 12-leads can be done (with the electrodes in the same place).