


Initial: 2/13/08
Reviewed/revised: 10/13/10
Revision: 1

**MILWAUKEE COUNTY EMS  
PRACTICAL SKILL  
ECG TRANSMISSION**

Approved by: Ronald Pirrallo, MD, MHSA
Signature:
Page 1 of 1

<b>Purpose:</b>	<b>Indications:</b>		
To transmit 4- lead and 12-lead electrocardiograms	Any patient who requires cardiac monitoring		
<b>Advantages:</b>	<b>Disadvantages:</b>	<b>Complications:</b>	<b>Contraindications:</b>
Transmits ECG to medical control and enables faxing 12-lead ECG to receiving hospital	None	None	None



**NOTES:**

- The 12-lead will remain in the Rosetta storage as long as the Rosetta device is powered on or for 20 minutes after the device turns itself off
- Paramedics can move the Rosetta and UHF radio to transmit from an area with better reception after disconnecting the Rosetta from the Zoll