


Initial: 3/1/15
Reviewed/revised:
Revision:

**MILWAUKEE COUNTY EMS
PRACTICAL SKILL
TEMPERATURE MEASUREMENT**

Approved: M. Riccardo Colella, DO, MPH, FACEP
Page 1 of 1

Purpose: To measure patient's temporal artery temperature		Indications:	
Advantages: Contact free measurement Digital and/or audible instant readout More precise than other methods Better tolerated than other methods	Disadvantages: None	Complications: None	Contraindications: None



NOTES:

- Thermometer may not be as accurate if used outside the ambient temperature range of 10°C/ 50°F and 40°C/ 104°F or if relative humidity is higher than 85%